

Surgical Guide

ennoia^{vet}

Personalized Engineering for Surgery



USE OF THE DEVICE

- Cutting Guide
- Reduction Guide
- Drilling Guide
- Resection Guide



MATERIAL

Polymer or titanium




STERILISATION

Autoclavable device



COMPLETION TIME

2 to 4 weeks



ennoia ensures for you the production of all types of surgical guides from DICOM data and surgical planning validated by you.

Cutting, reduction, drilling or resection guides, suitable for surgical tools (saw blade, drill bit, etc.).

Single-use product for a single patient.